BEZJEHLOVÝ VSTUP

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GB

NEEDLELESS PORT

= Sterilised using ethylene oxide	= Do Not Reuse	= Batch Code
= Non-pyrogenic	= Read instructions for use	= Use by
= Do not Resterilize	= Manufacturer	= Do not use if package is damaged
= Catalogue Number	= No natural latex included	= Does not contain di (2-ethylhexyl)-phthalate (DEHP)

Materials used: PC BPA-free / lipid resistant, Silicone

Compatibility: With medical devices equipped with standard Luer connectors in conformity with standard ČSN EN 1707 or ČSN EN 20594-1.

Contraindication: Allergic reactions to any material used for production of the product concerned.

Adverse effects: None known.

Warning: Do not use the needleless port if the packing is opened or damaged. Do not use the needleless port after the expiry date. Repeated sterilisation may cause increasing fragility or formation cracks.

Storage: The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

Caution:

- Replacement of the needleless port shall be performed always after 7 days or 100 connections at maximum.
- In order to prevent damage to the needleless port, use only the standard connection Luer which is in accordance with harmonized standards (non-standard syringes, needles, closures or connectors will cause damage to the needleless port). Do not angle.
- Usable life of the needleless port is 5 years from the date of production.

Intended use:

The needleless port is a closed safety system for use without a needle, ensuring microbiological safety, which is designed for repeated application as well as aspiration from intravenous accessories (syringes, infusion sets, cannulas, taps, catheters, extension hoses) with Luer type connection, in infusion therapy, injection therapy, sample taking, chemotherapy and it can also be used in magnetic resonance examinations.

Priming volume: 0,09 ml

Flow Rate Averages: 1 psi: 550 ml/min = gravitational flow
3 psi: 1000 ml/min
5 psi: 1300 ml/min

Instructions for use:

- Under aseptic conditions, remove the sterile needleless port from the packing.
- Prior to application, disinfect the silicone membrane surface by wiping with alkohol disinfection.
- Wait at least 30 seconds until the disinfection solution has dried out.
- Carefully connect (a common luer syringe or a luer lock connector) by turning clockwise to ensure a safe connection.
- Disconnecting the needleless port is carried out by turning anticlockwise.

NOTE:

It is not necessary to use gloves while using the needleless port.

After use, products should be disposed of in line with prescribed procedures for disposal of medical (infectious) waste in accordance with local regulations.

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