8. listopadu 2018

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November 08, 2018

At Halyard, we seek to constantly improve upon the quality of healthcare products and services available in the marketplace. We strive to deliver clinically-superior products with remarkable service to improve the wellbeing of the people we touch every day. In accordance with this vision, we want to inform our customers of additional chemical and chemotherapy drug testing performed on our **STERLING\* Nitrile Powder-Free Exam** Gloves. In the US, the Food and Drug Administration requires submission and clearance of a premarket notification for medical exam gloves that are intended for use with chemotherapy drugs and requires glove products to be labeled with drugs tested and breakthrough times. Package labeling for Halyard exam gloves sold in the US meet this labeling requirement.

Although additional chemotherapy drug and chemical testing is not required, Halyard believes the more informed our customers are, the better they will be able to understand the appropriate use of our products. The following chemotherapy drugs and chemicals have been tested according to EN374-1 (where applicable), ASTM D6978 (chemotherapy drugs) and ASTM F739 (chemicals). This data is provided for your information but is not intended to expand the intended use of the products beyond the scope of what has been cleared by the FDA or covered by the safety and performance requirements of the CE certification and European Norms.

If you have any questions, please do not hesitate to contact your Halyard qualified sales representative for more information.

Halyard gloves are tested by accredited independent laboratories: Akron Rubber & Development Lab (US), Centexbel (Belgium)  
\*Registered Trademarks or Trademarks of O&M Halyard or its affiliates

Copy 01382

*Disclaimer: Document to be used exclusively within EMEA*

**Chemotherapy Testing Information - Permeation testing per ASTM D6978**

<b>Chemotherapy Drugs Tested</b>	<b>Breakthrough time (minutes) Sterling* Nitrile</b>
Arsenic Trioxide (1 mg/ml)	>240
Azacitidine (25 mg/ml)	>240
Bendamustine (5 mg/ml)	>240
Bortezomib (1 mg/ml)	>240
Bleomycin sulfate (15 mg/ml)	>240
Busulfan (6 mg/ml)	>240
Carboplatin (10 mg/ml)	>240
Carfilzomib (10 mg/ml)	>240
Carmustine (3.3 mg/ml)	<b>do not use</b>
Cetuximab (Erbix) (2 mg/ml)	>240
Cyclophosphamide (20 mg/ml)	>240
Cytarabine HCL (100 mg/ml)	>240
Dacarbazine (10 mg/ml)	>240
Decitabine (5 mg/ml)	>240
Docetaxel (10 mg/ml)	>240
Doxorubicin HCL (2 mg/ml)	>240
Ellence (2 mg/ml)	>240
Eribulin Mesylate (0.5 mg/ml)	>240
Etoposide (20 mg/ml)	>240
Fludarabine (25 mg/ml)	>240
Fluorouracil (50 mg/ml)	>240
Fulvestrant (50 mg/ml)	>240
Gemcitabine (38 mg/ml)	>240
Idarubicin (1 mg/ml)	>240
Ifosfamide (50 mg/ml)	>240
Irinotecan (20 mg/ml)	>240
Mechlorethamine HCL (1 mg/ml)	>240
Melphalan (5 mg/ml)	>240
Methotrexate (25 mg/ml)	>240
Mitomycin (0.5 mg/ml)	>240
Oxaliplatin (2 mg/ml)	>240
Paraplatin (10 mg/ml)	>240
Pemetrexed (25 mg/ml)	>240
Pertuzumab (30 mg/ml)	>240
Rituximab (10 mg/ml)	>240
Temsirolimus (25 mg/ml)	>240
ThioTEPA (10 mg/ml)	<b>do not use</b>
Topotecan HCL (1 mg/ml)	>240
Trastuzumab (21 mg/ml)	>240

Halyard gloves are tested by accredited independent laboratories: Akron Rubber & Development Lab (US), Centexbel (Belgium)

\*Registered Trademarks or Trademarks of O&M Halyard or its affiliates

Copy 01382

*Disclaimer: Document to be used exclusively within EMEA*

Chemotherapy Drugs Tested	Breakthrough time (minutes) Sterling* Nitrile
Trisenox (1 mg/ml)	>240
Vinblastine (1 mg/ml)	>240
Vincristine Sulfate (1 mg/ml)	>240
Vinorelbine (10 mg/ml)	>240
Zoledronic Acid (0.8 mg/ml)	>240

**Chemical Testing Information- ASTM F739/EN374-1(where applicable)**

Chemicals Tested	CAS #	Breakthrough time (minutes) Sterling* Nitrile
Adhesive Tape Remover Pads		37.3
Alcoholic Bouin's Fixative		42.0
Alcoholic Eosin		8.0
Bleach	7681-9	> 480
Bouin's Fixative		> 480
B-Plus Fixative		> 480
CarNoy's Fluid		0
Cidex 14-Day		> 480
Cidex OPA		86.7
Decalcifier I		> 480
EM Fixative		> 480
Eosin Stain		42.0
<b>Ethidium Bromide (1%)</b>		<b>&gt;480</b>
Ethyl Alcohol (95%)	64-17-5	8.0
Formalin Alcohol (tested for Formaldehyde)		20.0
Formalin Alcohol (tested for Reagent Alcohol)		16.0
Formalin, 10% Buffered (tested for Formaldehyde)	50-00-0	> 480
Formalin, 10% Buffered (tested for Methanol)	50-00-0	> 480
Giemsa Stain Solution		19.3
Haematoxylin Stain (Harris Solution)		180.0
<b>Hydrochloric Acid, 30%</b>		<b>&gt; 480</b>
Hydrochloric Acid, 35.5%	7647-01-0	91.0
Hydrogen Peroxide, 3%	7722-84-1	> 480
Isopropyl Alcohol, 70%	67-63-0	30.9
Methanol	67-56-1	8.0
Methyl Methacrylate (Bone Cement)	80-62-6	0
Monse's Solution	1310-45-8	> 480
OxyCide Concentrate		40.0
OxyCide Ready to Use		> 480

Halyard gloves are tested by accredited independent laboratories: Akron Rubber & Development Lab (US), Centexbel (Belgium)  
 \*Registered Trademarks or Trademarks of O&M Halyard or its affiliates

Copy 01382

*Disclaimer: Document to be used exclusively within EMEA*

<b>Chemicals Tested</b>	<b>CAS #</b>	<b>Breakthrough time (minutes) Sterling* Nitrile</b>
Permaslip Mounting Medium & Liquid Cover Slip		0
Rapid Bone Decalcifier		> 480
Reagent Alcohol	64-17-5	10.7
Resert XL HLD		13.3
<b>Sodium Hydroxide (40%)</b>	1310-73-2	> 480
<b>Sulfuric Acid (50%)</b>	7664-93-9	> 480
Trichloroacetic Acid, (10%)	76-03-9	> 480
Xylene, (99%)	1330-20-7	0
<b>Acetone</b>	<b>67-64-1</b>	<b>&lt; 5</b>
<b>Ethyl Acetate</b>	<b>141-78-6</b>	<b>&lt; 5</b>
<b>n-Heptane</b>	<b>142-85-5</b>	<b>&gt; 480</b>
<b>Sulphuric Acid (96%)</b>	<b>7664-93-9</b>	<b>&lt; 5</b>
<b>Nitric acid (65%)</b>	<b>7697-37-2</b>	<b>&lt; 5</b>
<b>Acetic acid (99%)</b>	<b>64-19-7</b>	<b>&lt; 5</b>
<b>Ammonium Hydroxide (25%)</b>	<b>1336-21-6</b>	<b>7</b>
<b>Hydrogen Peroxide (30%)</b>	<b>7722-84-1</b>	<b>24</b>
<b>Formaldehyde (37%)</b>	<b>50-00-0</b>	<b>&gt; 480</b>

Halyard gloves are tested by accredited independent laboratories: Akron Rubber & Development Lab (US), Centexbel (Belgium)  
 \*Registered Trademarks or Trademarks of O&M Halyard or its affiliates

Copy 01382

*Disclaimer: Document to be used exclusively within EMEA*