- GB SILICONE COATED LATEX URINARY FOLEY CATHETERS – SOFT LATEX
- FR SONDES URINAIRES DE FOLEY EN LATEX SOUPLE AVEC REVÊTEMENT SILICONE
- DE MIT SILIKON BESCHICHTETE FOLEY-BLASENKATHETER AUS LATEX – WEICHER LATEX
- IT CATETERI URINARI IN LATTICE FOLEY CON RIVESTIMENTO IN SILICONE - LATTICE MORBIDO
- ES SONDAS FOLEY URINARIAS DE LÁTEX CON REVESTIMIENTO DE LÁTEX FLEXIBLE
- PT CATETERES URINÁRIOS DE FOLEY EM LÁTEX REVESTIDOS A SILICONE – LÁTEX MACIO
- NL LATEX FOLEY-URINEKATHETERS MET SILICONENCOATING – ZACHTE LATEX
- GR ΟΥΡΗΤΗΡΙΚΟΙ ΚΑΘΗΤΗΡΕΣ FOLEY ΛΑΤΕΞ ΜΕ ΕΠΙΣΤΡΩΣΗ ΣΙΛΙΚΟΝΗΣ – ΜΑΛΑΚΟ ΛΑΤΕΞ
- RU КАТЕТЕРЫ ФОЛЕЯ ИЗ ЛАТЕКСА С СИЛИКОНОВЫМ ПОКРЫТИЕМ ДЛЯ КАТЕТЕРИЗАЦИИ МОЧЕВОГО ПУЗЫРЯ — МЯГКИЙ ЛАТЕКС
- PL CEWNIKI MOCZOWE FOLEYA Z LATEKSU POKRYWANE SILIKONEM — MIĘKKI LATEKS
- CZ SILIKONEM POTAŽENÝ LATEXOVÝ MOČOVÝ FOLEY KATÉTR – MĚKKÝ LATEX
- TR SİLİKON KAPLAMA LATEKS ÜRİNER FOLEY KATETERİ - YUMUŞAK LATEKS
- FI PEHMEÄSTÄ LATEKSISTA VALMISTETUT SILIKONIPÄÄLLYSTETYT FOLEY-VIRTSAKATETRIT
- SE URINKATETER AV FOLEYTYP I SILIKONBELAGD LATEX – MJUK LATEX
- DK SILIKONEBELAGTE FOLEY-URINKATETRE AF LATEX – BLØD LATEX
- NO FOLEY-URINKATETERE I LATEKS MED SILIKONBELEGG – MYK LATEKS
- CN 硅涂层乳胶 FOLEY 导尿管 — 软乳胶
- BG ЛАТЕКСОВИ УРИНАРНИ ФОЛЕЙ КАТЕТРИ СЪС СИЛИКОНОВО ПОКРИТИЕ – МЕК ЛАТЕКС
- HU SZILIKON BEVONATÚ LATEX HÚGYÚTI FOLEY-KATÉTEREK – PUHA LATEX
- LT SILIKONU DENGTI LATEKSINIAI „FOLEY“ ŠLAPIMO KATETERIAI – IŠ MINKŠTO LATEKSO
- RO CATETERE URINARE FOLEY DIN LATEX SILICONAT – LATEX MOALE
- RS URINARNI FOLEY KATETERI OD LATEKSA OBLOŽENI SILIKONOM – MEKANI LATEKS
- العربية قساطر فولى البولوية من اللاتكس المغطى بالسيليكون- اللاتيكس الناعم

SILICONE COATED LATEX URINARY FOLEY CATHETERS – SOFT LATEX

Coloplast A/S
Holtedam 1
3050 Humlebaek – Denmark

REF.: SH2203
Date of issuance: September 2017
Revision Number: 2

SILICONE COATED LATEX URINARY FOLEY CATHETERS – SOFT LATEX

For Rx Only

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Please read all instructions before using this device.

DESCRIPTION

Silicone coated latex urinary catheter.

See the pack label for the length, diameter, balloon volume, material, method of sterilization and other specific characteristics.

The choice of the size, balloon volume and material of catheter is under the responsibility of the healthcare professional based on the patient status and planned indwelling time.

INDICATIONS

- Urethral urinary catheterization.

CONTRAINDICATIONS

Same as for urethral urinary catheterization.

Caution: *these products contain natural latex, which may cause allergic reactions.*

The evaluation of the allergic background of a patient is the health care professional's responsibility.

DURATION OF USE

Warning: These Soft latex Foley catheters cannot be left in place more than 7 days and must not be replaced immediately by the same or an identical catheter.

INSTRUCTIONS FOR USE

1) Urethral Insertion

- Choose a catheter of the appropriate size.
- Lubricate the catheter using a water-based gel and insert it according to the normal urethral catheterization technique, observing the usual procedures of asepsis.

N.B: Prior to catheterization, it is usual practice to check that the valve and balloon are functioning properly by inflating and then deflating the balloon.

Warning: *Paediatric catheters of diameter 08 and 10 FR/CH include a stylet that facilitates their insertion. They therefore have a closed end.*

Before insertion: Make sure that the stylet can move in the catheter and examine the end of the catheter to ensure that the stylet is positioned perfectly inside the catheter and does not come out of an eye.

- Make sure the catheter is positioned correctly by checking that there is urine at the external connector.

Warning: *After insertion of the paediatric catheters, withdraw the stylet using the orange grip, if necessary.*

2) Inflating the balloon

- Inflate the balloon with sterile water to the volume indicated on the package label. Connect a syringe without a needle to the anti-reflux valve and inject a volume of solution adapted to the volume of the balloon. Quickly disconnect the syringe once the balloon has been inflated and the catheter immobilized. Check that the balloon is inflated by pulling the catheter body gently – the catheter should offer some resistance.
- Connect the catheter to a urine bag.
- Check there is no leak at the connection point, and that the urine flows freely in the tube.

WARNINGS/PRECAUTIONS:

- This type of device must only be used by trained and experienced professionals.
- Any use other than stated indications is under the responsibility of the physician.
- If the catheter needs to be secured, the adhesive must be applied to the connector.
- Do not use petroleum-based lubricants with latex catheters.
- Do not inflate the balloon beyond the maximum value indicated.
- Do not clamp the catheter. Use a plug if necessary.
- Reuse of this single use product may create a potential risk to the user. Reprocessing, cleaning, disinfection and sterilization may compromise product characteristics which in turn create an additional risk of physical harm to or infection of the patient.

FOR THE PAEDIATRIC CATHETERS DIAMETER 08 AND 10 FR/CH:

- Do not inflate the balloon to more than 3 ml.
- Withdraw the stylet carefully so as not to move the catheter.

ADVERSE EVENTS

Several adverse events have been described with the use of balloon catheters.

Some are related to the patient's conditions, the others to the procedure or the device:

- Related to the patient: bladder irritation symptoms, pain, urinary tract infection, incrustation and stone formation,
- Related to the procedure: urinary tract trauma
- Related to the device: leakage, balloon burst or deflation.

FOLLOW-UP

Regular monitoring should be implemented to ensure that no side effect occurs, that the catheter is working properly and in particular that it is draining properly.

Monitoring should also check that the balloon is still inflated and the catheter correctly placed.

ADVICE TO THE PATIENT

- Patients should be educated on their indwelling catheter and the need for a regular monitoring.
- They should be advised to inform the attending physician immediately if any anomaly or dysfunction is noted.
- Ensure scrupulous local hygiene, cleaning the genital and anal areas every morning and evening with soap and water.

REMOVAL

This is carried out by simple traction after deflating the balloon using a syringe connected to the valve of the inflation lumen.

Balloon folds can be prevented by deflating gently and progressively the balloon. If the patient feels pain when the catheter is removed, the balloon can be slightly re-inflated (make sure the balloon is correctly placed inside the bladder before re-inflating it) and deflated gently once again to remove the folds.

Warning: *In exceptional circumstances, it can be difficult or even impossible to deflate the balloon. In the event of this please refer to your local policies to resolve the situation.*

The physician is sole judge of the duration of placement of a catheter depending of the type of catheter, its indication and the patient medical condition.

Soft latex Foley catheters must not be left in place for more than 7 days. After use, if treatment is still needed, the catheter must not be replaced by the same or another Soft latex Foley catheter.

STORE AWAY FROM LIGHT IN A COOL AND DRY PLACE.

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Hongkong – Keymax Technology, Unit 2103, 21/F, Kodak House II, 39 Healthy Street East, North Point- Hong Kong
Phone :+852 31508913 / Fax :+852 31508412 – keymaxtechnology@yahoo.com.hk

Colombia – Ordenes y Suministros , Carrera 13 A No78-88, barrio El Retiro. Bogota
Fax : 57 16347244 / E-mail: ivonnesalcedob@hotmail.com

Cyprus – Cyprus Pharmaceutical, P.O. Box 1005, 11 King Paul Street, 1500 Nicosia
Phone : +357 2863 100 / Fax: +357 2675 136 – cpomeditec@papellinas.com

Czech Republic – Promedica Praha Grp , Novodvorska 136 – 142 00 Praha 4
Phone : +420 /221 595 111 / Fax +420 224 257 838 rozsafny@promedica-praha.cz

Dominican Republic – Santos & Ortiz Group, Calle Dr Angel Severo Cabral69, Ensanche Julieta- Santo Domingo
Phone : +1 809 979 92 44 – rufinosantosh@gmail.com

Egypt – Egyptian Group for Import and Export, 11 Ebad El Rahman St, behind Tourism Institute, Sheraton Heliopolis 11799 Cairo
Phone : +20 2268 7712 / Fax: +20 2268 7714 – shymaa.abdelrahman@egggroup-eg.com

Egypt – Egypt care for Supplies & Commercial Agencies
4H Dar El-Defaa El Gawy st., El Nozha st. takseem Asmaa Fahmy, Ard El Golf, Heliopolis flat10 Cairo .
Phone : +2 02 24 191471 / Fax : +2 02 24 153079 contact maimostafa@egypt-care.com

Estonia – AB Medical Group Estonia, Reti tee 12, Peetri kula, 75312, Rae vald, Harjumaa
Phone : +372 655 2310/610 6461 / Fax: +372 659 3250 – grete@abmedical.ee

Finland – MEDIPLAST FENNO OY, PO Box 153- FI- 01531 Vantaa
Phone : +358 9 276 360 / Fax: +358 9 2763 6199 – kari.bergman@fennomedical.fi

Greece – Anastasios Mavrogenis SA, Passov 15, GR11141 Athens
Phone : +30 2101010132 / Fax: +30 210 2020775 – m.mavrogeni@coloplast.gr

Hungary – C.PP Budapest Kft, Amerikai UT 33., H-1145 Budapest
Phone : +3614223700 / Fax: +36 1 422 37 08, – gergely.jancso@cppbudapest.hu

Iran – Bartar Co, Bldg N°33 Sarparast ST Cross, Bozorgmehr St, Sth Felestin, Tehran
Phone : + 98 21664 87360 – afn@group7.org

Israel – Sela Medical, 6 Ravnitzky St. Petach, 42799 Tikva
Phone : +972 3 9043433 / Fax: +972 3 9045474 – yuval@selamedical.co.il

Jordan – Al Eiman Drug store, PO Box 402 - 7th Circle, Airport Street Beside City Terminal, 11810 Amman
Phone : +962 6 581 0108 / Fax: +962 6 581 8850 – Samerissa@eiman-eds.com

Korea – Suntek Medical Systems, 3rd Fl. Zero BD- 84 Gil, 14 Teheran-Ro , Kangnam-Ku- Seoul 135-845
Phone : +82 2 5143214 / Fax +82 2 547 2446 – Park Chong Young – cypark99@chol.com

Kuwait - Ashraf & Co., P.O. Box 3555, 13036 Safat
Phone : +965 531 2960 / Fax: +965 533 6372 - consumamedical@ashraf.com

Lebanon - Saramed SARL, 1022 Ryad El Solh, 1107 2070- Beirut
Phone : +961 1577 111 / Fax: +961 1575 600 - m.hasswany@saramed.com

Lithuania - Skirgesa Ltd, Ekskavatorininku str. 1b, LT-52461 Kaunas
Phone : +370 37 45 77 46, 45 21 68 / +370 37 45 81 61 - saulius@skirgesa.lt

Malaysia - Pro Element Technologies , 18-34-A3 Persiaran Gurney, Gurney Tower-10250 Penang
Phone : +604 226 32 86 / Fax +604 227 3286 - info@pro-element.com.my

Malta - Pharma-Cos Ltd., Pharma-Cos House Triq. C. Portanier,
1732 Santa Venera HMR 11
Phone : +356 2144 1977 / Fax: +356 2144 1388 - mmifs@pharma-cos.com

Mexico - Geoworks Durango 225-204, Col. Roma - Mexico, D.F 06700
Phone : +52 55 5514 0381 - Alberto Favilla - geoworks@geoworks.com.mx

Morocco - Sahmed, 2, rue des Acacias & Boulevard Y El Mansour El Maarif, 20100 Casablanca
Phone : +212 2298 1177 / Fax: +212 2298 1178 - logistics@claesmedical.com

Norway - Endomed A.S., Prof. Birkenlands vei 34B - 1081 Oslo
Phone : +47 23 17 24 00 / Fax: +47 23 17 24 01 - endomed@endomed.no

Oman - Etihad Medical Services LLC., Way 6310, Bldg 868 Office 23- Walja, Ruwi-Muscat
Phone : +968 244 99087 - binualex@eigllcoman.com

Pakistan - F.W. Distributors, FW Building, Opposite Poonch House, Adamjee Road, Saddar - Rawalpindi
Phone : +92 515516729 / Fax: +92 515120061 - fw@fwdpak.com.pk

Peru - Famedical - AV. Aviacion 3143, Ofic.101, San Borja-Lima
Phone : +511 225 88 02 - importaciones@famedical.com.pe

Pologne - Skamex, Spolka z ogr. odp. Sp. k., Czestochowska 38/52, 93-121 Łódź
Phone : +48 42 677 14 11 / Fax: +48 42 672 40 10 -
Katarzyna.kowalczyk@skamex.com.pl

Qatar - Al Danah Medical Company W.L.L., Gate N°4 Naser Bin Khaled Complex, Salwa Rd, PO Box 14485
Doha - Phone : +974 4469 1122 - m.affi@aldanahmedical.com

Romania - Steril Romania, 3-5, Metalurgiei Blvd., Sector 4 041831, Bucharest
Phone : +4021 460 07 79 / Fax: +4021 460 04 52 - import@steril.ro

Russia - MPA Medical Partners, 8th Martha Str, 1 bld.12- Moscow 127083
Phone : +7-495 921 30 88 / info@mpamed.ru

Russia - UROMED, Vvdenskogo str., 1 bld. 1 - Moscow 117342
Phone : +7-495 783 68 11 / info@uromed-m.ru

Saudi Arabia - JAT - 80 Amr Ibn Al-Aas street, Al Wisham Area Office 308 , Riyadh
Phone : +966 114040470 / maman@janatarab.com

Senegal - Kima Health Partner , 89 Rue 03 Mamelles Elevage Ouakam - Dakar
Phone : +221 33 821 42 37 - kimakhp@gmail.com

Serbia - Magna Pharmacia, Resavska 60, 11000 Belgrade,
Phone : +381 11 32 84 744 - dejana.boskovic@magnapharmacia.rs

Slovakia - SEQ-MED S.R.O, Biskupicka 21 , 986 01 Filakova
Phone : +421.2 555 766 05 - seqr@stonline.sk

Slovenia - Apollonia d.o.o., Zaloska 155, 1000 Ljubljana
Phone : +386 1/548 35 10 / Fax: +386 1 540 03 97 - apollonia@amis.net

South Africa - Nu Angle, Postnet Suite 93, Private Bag X33, Rivonia
Phone : +27 114674145 - dawn@nuangle.co.za

Tunisia - Medico Chirurgical, 31 BD Hedi Saidi, 1005 Tunis
Phone : +216 71899 750 / Fax: +216 7189 7050 - medico.chirurgi@gnet.tn

Turkey - Intra Medikal Gerecler Pazarlama Ltd., Koresehitleri cad. No.48/50
34394 Istanbul
Phone : +90 212 274 4752 / Fax: +90 212 288 2082 - porges@intramedical.com.tr

United Arab-Emirates - Pharmatrade, P.O. Box 11397, Dubai
Phone : +971 4268 3339 / Fax: +971 4262 5735 - amgad@pharmatradeuae.com

Uruguay - H. Strattner - Neomed SRL ,Francisco Llambi 1556, Montevideo, Uruguay,
Phone : +598 2 622 8326/ Fax: +598 2 622 5567 - rocioabella@neomed.com.uy

USA - Coloplast Corp. Minneapolis, MN 55411. Phone: +1 612-337-7800 -
www.coloplast.com
For customer service, please call 1-800-533-0464 in USA; Outside USA, contact your local Coloplast representative

Vietnam - Toan Anh CO. LTD, Appt. 300, Room 218, District 1, Hôchiminh City
Phone : +848 39574066, ivf_vietnam@yahoo.com

Rx ONLYCaution:
Federal (USA) law restricts this device to sale by or on the order of a physician**STERILE****R**

Sterilized using irradiation

Do not re-use

Use By

Date of manufacture

Caution, Consult accompanying documents

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