<b>GB</b>	BALLOON URINARY CATHETERS
<b>FR</b>	SONDES URINAIRES À BALLONNET
<b>DE</b>	BALLONBLASENKATHETER
<b>IT</b>	CATETERI URINARI A PALLONCINO
<b>ES</b>	SONDAS VESICALES CON BALÓN
<b>PT</b>	CATETERES URINÁRIOS DE BALÃO
<b>NL</b>	URINAIRE CATHETERS MET BALLON
<b>GR</b>	ΟΥΡΗΘΡΙΚΟΙ ΚΑΘΗΤΗΡΕΣ ΜΕ ΜΠΑΛΟΝΙ
<b>RU</b>	БАЛЛОННЫЕ МОЧЕВЫЕ КАТЕТЕРЫ
<b>PL</b>	CEWNIKI MOCZOWE Z BALONEM
<b>CZ</b>	BALÓNKOVÉ KATÉTRY PRO DRENÁŽ MOČOVÉHO MĚCHÝŘE
<b>TR</b>	BALON ÜRİNER KATETERLER
<b>FI</b>	VIRTSAPUTKEN PALLOKATETRIT
<b>SE</b>	BALLONGURINKATETRAR
<b>DK</b>	URINVEJSKATETER MED BALLON
<b>NO</b>	URINVEISKATETRE MED BALLONG
<b>CN</b>	球囊导尿管 产品描述
<b>BG</b>	БАЛОННИ УРИНАРНИ КАТЕТРИ
<b>HU</b>	HÚGYÚTI BALLONKATÉTEREK
<b>LT</b>	BALIONINIAI ŠLAPIMO KATETERIAI
<b>RO</b>	CATETERE URINARE CU BALON
<b>RS</b>	URINARNI BALON KATETERI
العربية	تعليمات الاستخدام

## BALLOON URINARY CATHETERS

0543

**Coloplast A/S**  
Holtedam 1  
3050 Humlebæk - Denmark

REF.: SH2115  
Date of issuance: March 2015  
Revision Number: 2

GB

INSTRUCTIONS FOR USE SH2115

## BALLOON URINARY CATHETERS

### DESCRIPTION

Silicone, latex, PVC or NEOPLEX® balloon urinary catheter (with or without hydrophilic coating: see label). See the pack label for the length, diameter, balloon volume, material, method of sterilization and other specific characteristics.

The choice of the size, balloon volume and material of catheter is under the responsibility of the healthcare professional based on the patient status and planned indwelling time.

### INDICATIONS

- Foley catheters: urethral urinary catheterization.
- Only straight 2-way Foley catheters with a maximum balloon volume of 15 ml may be used for the supra-pubic approach (except for ribbed catheters).
- 3-way Foley catheters: urethral urinary catheterization and postoperative bladder irrigation-lavage.
- Prostate catheters:
  - short-term drainage of bladder urine,
  - postoperative bladder irrigation-lavage,
  - after prostate surgery: haemostasis of the prostatic fossa.

### CONTRAINDICATIONS

Same as for urethral urinary catheterization and supra-pubic bladder drainage, and generally, known allergic reactions due to the device material (e.g. latex).

Where indicated, some products contain latex.

**Caution:** these products contain natural latex, which may cause allergic reactions.

The evaluation of the allergic background of a patient is the health care professional's responsibility.

### INSTRUCTIONS FOR USE

#### 1) Urethral Insertion

- Choose a catheter of the appropriate size.

Lubricate the catheter using a water-based gel and insert it according to the normal urethral catheterization technique, observing the usual procedures of asepsis.

**N.B:** for catheters with a hydrophilic coating: just before insertion, moisten the catheter for 15 seconds with sterile water or physiological saline solution to activate the lubrication properties of the hydrophilic coating. Ensure that this coating is kept moist during insertion. Do not use additional lubricants.

**N.B:** Prior to catheterization, it is usual practice to check that the valve and balloon are functioning properly by inflating and then deflating the balloon.

- Some catheters have an open end, which enables them to be inserted with a guidewire.

- Insert the catheter according to the usual procedure, by advancing it on the guidewire.

- After checking that the catheter is positioned properly, remove the guidewire.

**Warning:** paediatric catheters of diameter 06, 08 and 10 FR/CH may include a stylet that facilitates their insertion. They therefore have a closed end.

**Before insertion:** Make sure that the stylet can move in the catheter and examine the end of the catheter to ensure that the stylet is positioned perfectly inside the catheter and does not come out of an eye.

Paediatric catheters with guidewire:

- For catheters diameter 06 FR/CH, a diameter of 0.025 inch
- For catheters diameter 08 FR/CH, a diameter of 0.032 inch
- For catheters diameter 10 FR/CH, a diameter of 0.038 inch

- Make sure the catheter is positioned correctly by checking that there is urine at the external connector.

**Warning:** after insertion of the paediatric catheters, withdraw the stylet using the orange grip, if necessary.

#### 2) Supra-pubic insertion

- Follow the usual procedure for supra-pubic catheter insertion.

- Supra-pubic drainage must only be implemented when the bladder is full.

- Respect the rules of surgical asepsis.

#### 3) Inflating the balloon

- Inflate the balloon with sterile water to the volume indicated on the package label. Connect a syringe without a needle to the anti-reflux valve and inject a volume of solution adapted to the volume of the balloon. Quickly disconnect the syringe once the balloon has been inflated and the catheter immobilized. Check that the balloon is inflated by pulling the catheter body gently -the catheter should offer some resistance.

- Connect the catheter to a urine bag.

- Check there are no leaks at the connection point, and that the urine flows freely in the tube.

### WARNINGS/PRECAUTIONS:

- This type of device must only be used by trained and experienced professionals.

- Any use other than stated indications is under the responsibility of the physician.

- If the catheter needs to be secured, the adhesive must be applied to the connector.

- To lubricate catheters without a coating, a water-based lubricant is recommended.

- Do not use petroleum-based lubricants with latex catheters.

- Do not use silicone oil with silicone catheters.

- Using an iodine-based irrigation-lavage solution may damage some silicone catheters.

- Do not inflate the balloon beyond the maximum value indicated.

- Do not clamp the catheter. Use a plug if necessary.

- Reuse of this single use product may create a potential risk to the user. Reprocessing, cleaning, disinfection and sterilization may compromise product characteristics which in turn create an additional risk of physical harm to or infection of the patient.

### PRECAUTIONS FOR THE PAEDIATRIC CATHETERS DIAMETER 06, 08 AND 10 FR/CH.

*Catheters diameter 06 FR/CH: Do not inflate the balloon to more than 1.5 ml.*

*Catheters diameter 08 and 10 FR/CH: Do not inflate the balloon to more than 3 ml.*

Withdraw the stylet carefully so as not to move the catheter.

### ADVERSE EVENTS

Several adverse events have been described with the use of balloon catheters.

Some are related to the patient's conditions, the others to the procedure or the device:

- Related to the patient: bladder irritation symptoms, pain, urinary tract infection, incrustation and stone formation.

- Related to the procedure: urinary tract trauma.

- Related to the device: leakage, balloon burst or deflation.

Adverse events specifically related to the use of supra-pubic catheters have been described.

- Related to the patient: same as above and hematuria, any type of skin irritation.

- Related to the procedure: the potential adverse events are those observed with supra-pubic drainage, particularly if one forgets to check whether the bladder is full before puncturing it, among which peritoneal perforation with or without bowel perforation, misplacement/displacement.

- Related to the device: same as above and migration of the catheter, catheter knotting.

### FOLLOW-UP

Regular monitoring should be implemented to ensure that no side effect occurs, that the catheter is working properly and in particular that it is draining properly.

Monitoring should also check that the balloon is still inflated and the catheter correctly placed.

### ADVICE TO THE PATIENT

- Patients should be educated on their indwelling catheter and the need for a regular monitoring.

- They should be advised to inform the attending physician immediately if any anomaly or dysfunction is noted.

- Ensure scrupulous local hygiene, cleaning the genital and anal areas every morning and evening with soap and water.

### REMOVAL

This is carried out by simple traction after deflating the balloon using a syringe connected to the valve of the inflation lumen.

Especially for silicone catheters, balloon folds can be prevented by deflating gently and progressively the balloon. If the patient feels pain when the catheter is removed, the balloon can be slightly re-inflated (make sure the balloon is correctly placed inside the bladder before re-inflating it) and deflated gently once again to remove the folds.

**Warning:** in exceptional circumstances, it can be difficult or even impossible to deflate the balloon. In the event of this please refer to your local policies to resolve the situation.

The physician is sole judge of the duration of placement of a catheter depending of the type of catheter, its indication and the patient's medical condition.

FOLYSIL® silicone catheters must not be left in place for more than 30 days.

Foley NEOPLEX® and PVC catheters must not be left in place for more than 2 weeks.

Foley latex catheters must not be left in place for more than 7 days.

Prostatic catheters must not be left in place for more than 7 days.

STORE AWAY FROM LIGHT IN A COOL AND DRY PLACE.

DO NOT USE IF PACKAGE IS DAMAGED.

DO NOT RESTERILIZE.

DISCARD PRODUCT AFTER USE.

### For Rx Only

**Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.**

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**Pakistan** – F.W. Distributors, FW Building, Opposite Poonch House, Adamjee Road, Saddar - Rawalpindi  
Phone : +92 515516729 / Fax: +92 515120061 – [fw@fwdpak.com.pk](mailto:fw@fwdpak.com.pk)

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**USA** – Coloplast Corp. Minneapolis, MN 55411. Phone : +1 612-588-4685 - [www.coloplast.com](http://www.coloplast.com)  
For customer service, please call 1-800-533-0464 in USA;Outside USA, contact your local Coloplast representative

**Venezuela** – Althea Medical CA, Av. Romulo Gallegos, Torre Exagon Piso 9, No. 92,  
El Marques, Caracas  
Phone: +58 212 242 0234 / Fax: +58 212 242 0358 - [altheamedical@gmail.com](mailto:altheamedical@gmail.com)

**Vietnam** – Toan Anh CO. LTD, Appt. 300, Room 218, District 1, Hôchiminh City  
Phone: +848 39574066, [jvf\\_vietnam@yahoo.com](mailto:jvf_vietnam@yahoo.com)

**Rx ONLY** Caution:  
Federal (USA) law restricts this device to sale by or on the order of a physician

**STERILE** EO Sterilized using ethylene oxide

 Do not re-use

 Use By

 Date of manufacture

 Caution, Consult accompanying documents

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The Coloplast logo is a registered trademark of Coloplast A/S.

© [2015-03.]

All rights reserved. Coloplast A/S, 3050 Humlebaek, Denmark.

Signature Page for VV-0095337 v1.0

Approved by	FRJAC Jean Ann Caminade RMF Specialist & prod support Technical / Specialist 24-Jun-2015 07:32:17 GMT+0000
-------------	---

Signature Page for VV-0095337 v1.0

Document Owner: FRJAC Jean Ann Caminade