- GB** FOLYSIL® LT AND X-TRA LONG-TERM CATHETERISATION
- FR** FOLYSIL® LT et X-TRA
SONDAGE A DEMEURE LONGUE DURÉE
- DE** FOLYSIL® LT UND X-TRA
LANGZEIT-KATHETERISIERUNG
- IT** FOLYSIL® LT E X-TRA
CATETERISMO A LUNGO TERMINE
- ES** FOLYSIL® LT Y X-TRA
SONDAJE DE LARGA DURACIÓN
- PT** FOLYSIL® LT E X-TRA
CATETERIZAÇÃO PROLONGADA
- NL** FOLYSIL® LT EN X-TRA
VOOR LANGDURIGE CATHETERISATIE
- GR** FOLYSIL® LT KAI X-TRA
ΚΑΘΗΤΗΡΙΑΣΜΟΣ ΜΑΚΡΑΣ ΔΙΑΡΚΕΙΑΣ
- RU** FOLYSIL® LT И X-TRA КАТЕТЕР ДЛЯ
ДОЛГОВРЕМЕННОЙ УСТАНОВКИ
- PL** FOLYSIL® LT I X-TRA
DO DŁUGOTRWAŁEGO CEWNIKOWANIA
- CZ** FOLYSIL® LT A X-TRA
DLUHODOBÁ KATETRIZACE
- TR** FOLYSIL® LT VE X-TRA
UZUN SÜRELİ KATETER UYGULAMASI
- FI** FOLYSIL® LT JA X-TRA
PITKÄAIKAINEN KATETROINTI
- SE** FOLYSIL® LT OCH X-TRA
LÅNGTIDSKATETERISERING
- DK** FOLYSIL® LT OG X-TRA
KATETER Å DEMEURE
- NO** FOLYSIL® LT OG X-TRA
LANGTIDSKATETERISERING
- CN** FOLYSIL® LT 與 X-TRA 長期放置導尿管
- BG** FOLYSIL® LT И X-TRA
ПРОДЪЛЖИТЕЛНА КАТЕТЕРИЗАЦИЯ
- HU** HOSSZÚ TÁVÚ KATÉTEREZÉS FOLYSIL® LT
és X-TRA KATÉTERREL
- LT** „FOLYSIL® LT“ IR „X-TRA“ ILGALAIKĖS
KATETERIZACIJOS KATETERIAI
- RO** CATETERIZARE PE TERMEN LUNG -FOLYSIL® LT ŞI X-TRA
- RS** FOLYSIL® LT I X-TRA DUGOTRAJNA KATETERIZACIJA
- العربية** القساطر طويلة المدى X-TRA و FOLYSIL® LT

FOLYSIL® LT AND X-TRA LONG-TERM CATHETERISATION

0543

Coloplast A/S
Holtedam 1
3050 Humlebæk - Denmark

REF.: SH2082
Date of issuance: December 2014
Revision Number: 2

GB

INSTRUCTIONS FOR USE SH2082

FOLYSIL® LT AND X-TRA LONG-TERM CATHETERISATION

DESCRIPTION:

Silicone Foley catheter.

Following ref AA81, AA85 are supplied with a prefilled syringe (containing 10ml sterile water with 10% glycerine) for balloon inflation.

See the label on the box for the length, diameter, balloon volume and other specific characteristics.

INDICATIONS:

Urethral and supra-pubic bladder drainage.

CONTRAINDICATIONS:

- Same as for urethral urinary catheterization and supra-pubic bladder drainage.
- The evaluation of the allergic background of a patient is the health care professional's responsibility.

OPERATING PROCEDURE:

NOTE: Prior to using a balloon catheter, it is usual practice to check that the valve and balloon are functioning properly by inflating and then deflating the balloon.

1) Urethral Insertion

- Choose a catheter of the appropriate size.
- Lubricate the catheter using a water-based gel and insert it according to the normal urethral catheterization technique, observing the usual procedures of asepsis.
- Some catheters have an open end, which enables them to be inserted with a guidewire.
- Insert the catheter according to the usual procedure, by advancing it on the guidewire.
- After checking that the catheter is positioned properly, remove the guidewire.
- Make sure the catheter is positioned correctly by checking that urine appears at the external connector.

WARNING: the paediatric catheters of diameter 06, 08 and 10 CH/FR include a stylet that facilitates their insertion.

Before insertion: make sure that the stylet can move properly in the catheter and examine the end of the catheter to ensure that the stylet is positioned perfectly inside the catheter and does not come out of an eye.

After insertion, withdraw the stylet using the orange grip.

2) Supra-pubic insertion

- Follow the usual procedure for suprapubic catheter insertion.
- **Supra-pubic drainage must only be implemented when the bladder is full.**
- **Respect the rules of surgical asepsis.**

3) Balloon inflation

- If the catheter is supplied with a prefilled syringe (Folysil X-Tra), use the syringe to inflate the balloon to the volume indicated on the package label.
 - If there is no prefilled syringe, inflate the balloon with sterile water to the volume indicated on the package label. Do use a syringe without a needle.
 - Quickly disconnect the syringe once the balloon has been inflated and the catheter immobilized.
- Check that the balloon is inflated by pulling the catheter body gently. The catheter should offer some resistance.

- Connect the catheter to a urine bag.

Check there are no leaks at the connection point, and that the urine flows freely in the tube.

PRECAUTIONS:

- This type of device must only be used by trained and experienced professionals.
- Any use other than stated indications is under the responsibility of the physician.
- The choice of the size, balloon volume and material of catheter is under the responsibility of the healthcare professional based on the patient status and planned indwelling time.
- If the catheter needs to be secured, the adhesive must be applied to the connector.
- To lubricate catheters, a water-based lubricant is recommended.
- Do not use silicone oil with silicone catheters.
- Using a iodine-based irrigation-lavage solution may damage some silicone catheters.
- Do not inflate the balloon beyond the maximum value indicated.
- Do not clamp the catheter. Use a plug valve or spigot if necessary.
- The prefilled syringe supplied with catheter (if applicable) is to be used only to inflate catheter balloons, and not for any other procedure. Not for injection!
- Do not reuse any liquid extracted from the catheter balloon when changing the catheter to refill the balloon of the newly placed catheter, especially not the glycerine solution provided in the prefilled syringe.
- Reuse of this single use product may create a potential risk to the user. Reprocessing, cleaning, disinfection and sterilization may compromise product characteristics which in turn create an additional risk of physical harm to or infection of the patient.

PRECAUTIONS SPECIFIC TO THE PAEDIATRIC CATHETERS OF DIAMETER CH/FR 06, 08 AND 10:

- CH/FR 06 catheters: do not inflate the balloon to more than 1.5 ml
- CH/FR 08 and 10 catheters: do not inflate the balloon to more than 3 ml.
- Withdraw the stylet carefully to prevent moving the catheter.

ADVERSE EVENTS:

Several adverse events have been described with the use of balloon catheters.

Some are related to the patient conditions, the others to the procedure or the device:

- Related to the patient: bladder irritation symptoms, pain, urinary tract infection, incrustation and stone formation,
- Related to the procedure: urinary tract trauma
- Related to the device: leakage, balloon burst or deflation.

Adverse events specifically related to the use of suprapubic catheters have been described.

- Related to the patient: same as above and, hematuria, any type of skin irritation.
- Related to the procedure: the potential adverse events are those observed with suprapubic drainage, particularly if one forgets to check whether the bladder is full before puncturing it, among which peritoneal perforation with or without bowel perforation, misplacement/displacement.
- Related to the device: same as above and migration of the catheter, catheter knotting.

FOLLOW-UP:

Regular monitoring should be implemented to ensure that no side effect occurs, that the catheter is functioning properly and in particular that it is draining properly.

Monitoring should also check that the balloon is still inflated and the catheter correctly placed.

ADVICE TO THE PATIENT

- Patients should be educated on their indwelling catheter or their suprapubic catheter and the need for a regular monitoring.
- They should be advised to inform the attending physician immediately if any anomaly or dysfunction is noted.
- Ensure scrupulous local hygiene, cleansing the genital and anal areas every morning and evening with soap and water, wiping away from the body.

REMOVAL:

This is carried out by simple traction after deflating the balloon using a syringe connected to the valve of the inflation lumen. In order to avoid balloon folds, gentle and slow deflation is recommended to empty the balloon.

Warning: In exceptional circumstances, it can be difficult or even impossible to deflate the balloon. In this event please refer to your local policies to resolve the situation.

Balloon folds can be prevented by deflating gently and progressively the balloon. If the patient feels pain when the catheter is removed, the balloon can be slightly re-inflated (make sure the balloon is correctly placed inside the bladder before re-inflating it) and deflated progressively once again to remove the folds. FOLYSIL® LT and X-Tra silicone catheters can be used for up to 12 weeks.

STORE AWAY FROM LIGHT IN A COOL AND DRY PLACE.

DO NOT USE IF PACKAGE IS DAMAGED.

DO NOT RESTERILIZE.

DISCARD PRODUCT AFTER USE.

DO NOT REUSE ANY PRODUCT PROVIDED IN THIS PACKAGING.

For Rx Only

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

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Hongkong – Keymax Technology, Unit 2103, 21/F., Kodak House II, 39 Healthy Street East,
North Point- Hong Kong
Phone: +852 31508913 / Fax: +852 31508412 – keymaxtechnology@yahoo.com.hk

Colombia – Ordenes y Suministros, Carrera 13 A No78-88, barrio El Retiro. Bogota
Fax: +57 16347244 / E-mail: ivonnesalcedob@hotmail.com

Cuba – Medicuba, Maximo Gomez n° 1, Habana Vieja, Ciudad Habana,
Phone: +53 7 862 40 61/62/63 / Fax: +53 7 866 85 16 medicuba@medicuba.sld.cu

Cyprus – Cyprus Pharmaceutical, P.O. Box 1005, 11 King Paul Street, 1500 Nicosia
Phone: +357 2863 100 / Fax: +357 2675 136 – cpomeditec@papellinas.com

Czech Republic – Promedica Praha Grp., Novodvorska 136 – 142 00 Praha 4
Phone: +420 /221 595 111 / Fax +420 224 257 838 petrusek@promedica-praha.cz

Egypt – Egyptian Group, 11 Ebad El Rahman St, behind Tourism Institute,
Sheraton Heliopolis 11799 Cairo
Phone: +20 2268 7712 / Fax: +20 2268 7714 – laila@egggroup-eg.com

Egypt – Egypt care for Supplies & Commercial Agencies
4H Dar El-Defaa El Gawy st., El Nozha st. takseem Asmaa Fahmy, Ard El Golf,
Heliopolis flat10 Cairo
Phone: +2 02 24 191471 / Fax: +2 02 24 153079 contact Maysar@egypt-care.com

Estonia – AB Medical Group Estonia, Reti tee 12, Peetri kula, 75312, Rae vald, Harjumaa
Phone: +372 655 2310/610 6461 / Fax: +372 659 3250 – grete@abmedical.ee

Finland – Fenno Medical OY, PO Box 153- FI- 01531 Vantaa
Phone: +358 9 276 360 / Fax: +358 9 2763 6199 – kari.bergman@fennomedical.fi

Greece – Anastasios Mavrogenis SA, Passov 15, GR11141 Athens
Phone: +30 2101010132 / Fax: +30 210 2020775 – m.mavrogeni@coloplast.gr

Hungary – C.P.P Budapest Kft, Amerikai UT 33, H-1145 Budapest HU
Phone: +3614223700 / Fax: +36 1 422 37 08, gergely.jancso@cppbudapest.hu

Iran – Arman Tandorost Co., Nr. 23, 5th floor Etminan Bldg 480, Tehran
Phone: +98 2166 413870 / Fax: +98 2166 412885 – atc@mavara.com

Israel – Sela Medical, 6 Ravnitzky St. Petach, 42799 Tikva
Phone: +972 3 9043433 / Fax: +972 3 9045474 – itai@selamedical.co.il

Jordan – Al Eiman Drug store, PO Box 402 – 7th Circle, Airport Street Beside City Terminal,
11810 Amman
Phone: +962 6 581 0108 / Fax: +962 6 581 8850 – Samerissa@eiman-eds.com

Korea – Suntek Medical Systems, 413 Woolim e Biz Ctr, 280-21 Seongsudung 2 – Ga
Seongdong-Ku, Seoul 133-120
Phone: +82 2 5143214 / Fax +82 2 547 2446 – Park Chong Young cypark99@chol.com

Kuwait – Ashraf & Co., P.O. Box 3555, 13036 Safat
Phone: +965 531 2960 / Fax: +965 533 6372 – consumamedical@ashraf.com

Lebanon – Saramed SARL, 1022 Ryad El Solh, 1107 2070- Beyrut
Phone: +961 1577 111 / Fax: +961 1575 600 – m.hasswany@saramed.com

Libya – Kema Medical, Ras-hassan Eman, Elgharezi Street -Tripoli
Phone: +218 1335 1197 / Fax: +218 1335 1199 Info@kema-ly.com

Lithuania – Skirgesa Ltd, Ekskavatorininku str. 1b, LT-52461 Kaunas
Phone: +370 37 45 77 46, 45 21 68 / +370 37 45 81 61 – saulius@skirgesa.lt

Malaysia – Pro Element Technologies, 18-34-A3 Persiaran Gurney, Gurney Tower-10250 Penang
Phone: +604 226 32 86 / Fax +604 227 3286 – info@pro-element.com.my

Malta – Pharma-Cos Ltd., Pharma-Cos House Triq. C. Portanier, 1732 Santa Venera HMR 11
Phone: +356 2144 1977 / Fax: +356 2144 1388 – mmifs@pharma-cos.com

Mexico – Geoworks Durango 225-204, Col. Roma – Mexico, D.F 06700
Phone: +52 55 5514 0381 – Alberto Favilla – geoworks@geoworks.com.mx

Morocco – Sahmed, 2, rue des Acacias & Boulevard Y El Mansour El Maarif, 20100 Casablanca
Phone: +212 2298 1177 / Fax: +212 2298 1178 – logistics@claesmedical.com

Norway – Endomed A.S., Prof. Birkenlands vei 34B – 1081 Oslo
Phone: +47 23 17 24 00 / Fax: +47 23 17 24 01, endomed@endomed.no

Oman – Muscat Pharmacy and Stores LLC, Po Box 438, 113 Muscat,
Phone: +968 24814501 / Fax: +968 24815201 – mphq@mpmct.com,
mansoor@nationalpharmacylc.com

Pakistan – F.W. Distributors, FW Building, Opposite Poonch House, Adamjee Road,
Saddar – Rawalpindi
Phone: +92 515516729 / Fax: +92 515120061 – fw@fwdpak.com.pk

Peru – Famedical - Calle Los Halcones 410, Ofic. 401 URB. Jardin, Surquillo- Lima
Phone: +511 723 2615 importaciones@famedical.com.pe

Pologne – Skamex, Spolka z ogr. odp. Sp. k., Czeszochowska 38/52, 93-121 Łódź
Phone: +48 42 677 14 11 / Fax: +48 42 672 40 10 – Katarzyna.kowalczyk@skamex.com.pl

Qatar – Ali Bin Ali, P.O. Box 75 Al Jelaia str 37 Bin Omran-behind Hamad Hosp-Centre Doha
Phone: +9744867871 / Fax: +9744882585 – shadi.sleiman@alibinali.com

Romania – Steril Romania, 3-5, Metalurgiei Blvd., Sector 4 041831, Bucharest
Phone: +4021 460 07 79 / Fax: +4021 460 04 52 – import@steril.ro

Saudi Arabia – JAT – 80 Amr Ibn Al-Aas street, Al Wisham Area Office 308 , Riyadh
Phone: +966 14040470 / maman@janatarab.com

Senegal – Kima Health Partner , 89 Rue 03 Mamelles Elevage Ouakam – Dakar
Phone: +221 33 821 42 37 – kimakhp@gmail.com

Serbia – Magna Pharmacia, Resavska 60, 11000 Belgrade,
Phone: +381 11 32 84 744 - dejana.boskovic@magnapharmacia.rs

Slovakia – SEQ-MED S.R.O , Biskupicka 21 SK-936 01 Filakova
Phone: +421.47 438 25 60 – seqr@stonline.sk

Slovenia – Apollonia d.o.o., Zaloska 155, 1000 Ljubljana
Phone: +386 1/548 35 10 / Fax: +386 1 540 03 97 - apollonia@amis.net

South Africa – Nu Angle, Postnet Suite 93, Private Bag X33, Rivonia
Phone: +27 114674145 - dawn@nuangle.co.za

Syria – Zarzour Co, P.O. Box 72, Al Tal, Aldaheia N.H Way- Damascus
Phone: +963 11 59 48 000 / Fax: +963 11 59 47 000 - az@zarzourmn.com

Tunisia – Medico Chirurgical, 31 BD Hedi Saidi , 1005 Tunis
Phone: +216 71899 750 / Fax: +216 7189 7050 - medico.chirurgi@gnet.tn

Turkey – Intra Medikal Gereeder Pazarlama Ltd., Koresehitleri cad. No.48/50
34394 Istanbul
Phone: +90 212 274 4752 / Fax: +90 212 288 2082 - porges@intramedical.com.tr

United Arab-Emirates – Pharmatrade, P.O. Box 11397, Dubai
Phone: +971 4268 3339 / Fax: +971 4262 5735 - amgad@pharmatradeuae.com

Uruguay – H. Strattner – Neomed SRL ,Francisco Llambi 1556, Montevideo, Uruguay,
Phone: +598 2 622 8326/ Fax: +598 2 622 5567 - mpantazoglu@neomed.com.uy

USA – Coloplast Corp. Minneapolis, MN 55411. Phone : +1 612-588-4685 - www.coloplast.com
For customer service, please call 1-800-533-0464 in USA;Outside USA, contact your local Coloplast representative

Venezuela – Althea Medical CA, Av. Romulo Gallegos, Torre Exagon Piso 9, No. 92,
El Marques, Caracas
Phone: +58 212 242 0234 / Fax: +58 212 242 0358 - altheamedical@gmail.com

Vietnam – Toan Anh CO. LTD, Appt. 300, Room 218, District 1, Hôchiminh City
Phone: +848 39574066, jvf_vietnam@yahoo.com

**Rx Only**

Caution:

Federal (USA) law restricts this device to sale by or on the order of a physician

STERILE | **EO**

Sterilized using ethylene oxide

Do not re-use

Use By

Date of manufacture

Caution, Consult accompanying documents

LOT

Batch Code

Keep away from sunlight

Do not use if package is damaged
 Ne pas utiliser si l'emballage est endommagé
 Bei beschädigter Verpackung nicht verwenden
 Non utilizzare il prodotto se la confezione è danneggiata
 No utilizar si el embalaje está dañado
 Não utilize se a embalagem estiver danificada
 Niet gebruiken als de verpakking is beschadigd
 Μη χρησιμοποιείτε το προϊόν αν η συσκευασία έχει υποστεί ζημιά
 Не использовать, если упаковка повреждена
 Nie używać, jeśli opakowanie jest uszkodzone
 V prípade poškodeného balení nepoužívajte
 Ambalaj hasarlıysa ürünü kullanmayın

Ei saa käyttää, jos pakkaus on vahingoittunut
 Använd inte om förpackningen är skadad
 Må ikke anvendes hvis emballagen er beskadiget
 Skal ikke brukes hvis emballasjen er skadet
 包装破损不可使用
 Да не се използва, ако опаковката е повредена
 Ne használja, ha a csomagolás nyitva van vagy sérült.
 Negalima naudoti, jei pakuotė pažeista
 A nu se utiliza dacă ambalajul este deteriorat
 Ne koristiti ako je paklet oštećen
 منع الاستخدام في حال تلف العبوة

Not manufactured with Natural Rubber Latex
 N'est pas fabriqué avec du latex naturel
 Nicht mit Naturkautschuklatex hergestellt
 Non prodotto con lattice di gomma naturale
 Este producto no esta fabricado con goma de latex natural
 Nao fabricado com borracha de latex natural
 Niet gemaakt uit natuurlijke rubberlatex
 Δεν είναι κατασκευασμένο από φυσικό λάτεξ
 Произведен без применения натурального каучукового латекса
 Do produkcji nie użyto lateksu naturalnego
 Vyrobeno bez přírodní gumy
 Doğal kauçuk lateks ile uretilmemiştir

Valmistettu ilman luonnonkumilateksia
 Produkten är inte tillverkad av aturgummitatex
 Inholder ikke naturgummitatex
 Ikke produsert med naturgummitatex
 不由天然乳胶制成
 Не е произведено с латекс от естествен каучук
 Nem tartalmaz természetes gumilateket
 Sudetyje nėra natūralaus latekso
 Nu este fabricat cu latex din cauciuc natural
 Nije proizvedeno od prirodne lateks gume
 تسببه الطبيعي في المطاطي اللاتكس يدخل لا

The Coloplast logo is a registered trademark of Coloplast A/S.

© [2014-12.]

All rights reserved Coloplast A/S, 3050 Humlebæk, Denmark.