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8IN035 Rev.01 01-2008 Rev.10 09-2019

EN Foley Catheter with Temperature Sensor INTENDED USE

The Foley Catheter with Temperature Sensor is to be used for drainage of the urinary bladder and simultaneous monitoring of temperature. Foley catheterization is recommended when continuous bladder drainage is required in the absence of voluntary urination, after surgery or trauma of the pelvic organs, in the presence of obstruction or paralysis, and as a means to evaluate urine output in patients having fluid management problems. It is sterile and is individually packaged. It is of all silicone construction. The balloon inflation system is terminated with a Luer tip syringe actuated, self-sealing valve. A temperature sensor is located near the distal tip of the catheter. The catheter size and balloon volume are designated on the catheter itself.

The temperature is taken by the 400 series thermistor probe which has been inserted into the third channel up to the drainage eye of the catheter. The temperature range is 25°C to 45°C, with a tolerance of ±0.1°C. The data is sent to a monitor for processing via a connector and the re-usable adapter cable.

READ ALL INSTRUCTIONS CAREFULLY BEFORE USE

The Foley Catheter with Temperature Sensor must only be connected to a medical device registered in the EU (CE marked), FDA or any other official medical registration authority, and for which this catheter is an authorized accessory!
The monitoring system must comply with EN 60601-1.

ADVERSE REACTIONS

The following have been reported to be associated with the use of Foley Catheter and/or temperature sensors, especially during prolonged use: bladder spasms, electric burns due to aberrant electro-surgical current, encrustations with urinary salts, kidney or bladder stones, septicemia, urethral discharge, urethritis, urethrocutaneous fistulae and urinary tract infection or irritation.

PRECAUTIONS

- Follow the recommended balloon inflation volume printed on the product and use sterile water for inflation.
- Select the appropriate size catheter using expert clinical and medical judgment.
- Lubricate catheter with a suitable water soluble lubricant prior to insertion using accepted medical techniques.
- Deflate the balloon prior to removal of the catheter.
- Do not use a needle to puncture the catheter wall for aspiration. Leaks may occur.
- Do not intertwine the cables, especially the monitor cables with the electro-surgical unit’s cables.
- The operation of the patient temperature monitor may be temporarily affected during electro-surgical activations. Unusual temperature readings should be checked.

DIRECTIONS FOR USE

- Use aseptic technique in patient preparation and catheter care.
- Remove the Foley Catheter with Temperature Sensor, without stretching, from its sterile package.
- Visually inspect the catheter for any mechanical damage. If a trfial inflation of the balloon is carried out, ensure the balloon is empty prior to inflation. Inspect compatibility of sensor, cable and monitor system. Do not use if sensor or cable appear to be cracked, broken or twisted apart.
- Lubricate the catheter with a suitable water soluble lubricant and insert using acceptable medical techniques.
- Make sure that the catheter is draining properly before inserting the balloon.
- To inflate the balloon, insert a Luer tip syringe filled with the recommended amount of sterile water into the valve, i.e. 3cc or 5cc. Depress the plunger. The valve will seal when the syringe tip is withdrawn.
- Attach drainage tubing and bag.
- Secure the catheter, drainage tubing and bag in place following accepted medical technique.
- Connect the temperature probe to the appropriate reusable cable. Connect the cable to the patient monitor. Secure the cable using the drape clip.
- Follow the directions for use of the patient temperature monitor instrument.
- Prior to removing the catheter, deflate the balloon by inserting the Luer tip of an empty syringe into the valve and aspirating the solution completely.
- If balloon does not collapse after attempting to perform #11 above, sever the valve arm. If balloon still does not deflate, call physician.
- Remove the catheter following accepted medical techniques.
- Disconnect the probe (at the connector) and discard probe. To disconnect, grasp both connectors firmly and pull. Do not pull on cable or wire.

POTENTIAL COMPLICATIONS

If difficulty is encountered aspirating the balloon with a syringe, a rare and infrequently reported event, the channel of the catheter with the color coded valve should be cut with sharp scissors at the bifurcation, or the balloon should be ruptured according to established procedures reported in medical literature. Should it be necessary to rupture the balloon, care must be taken to remove all fragments from the patient.
Incorrectly positioned catheters can be cause urethral damage if the balloon is inflated within the urethra. Irritation of the urethral mucosa, electrical burns due to faulty electric currents or circuits, blockage of the catheter due to encrustation and catheter induced infections are documented complications. The patient should be routinely monitored in accordance with accepted procedures. The catheter should be removed after a suitable interval, can be longer than 30 days, as determined by a physician or other suitably qualified personnel.

CAUTION

- For urological use only. Federal law restricts this device to sale by or on the order of a physician (USA & Canada).
- Use only a water based lubricant.
- Always inflate the balloon with sterile water.
- Never clamp the catheter. A catheter plug or valve may be used if required.
- Never use a syringe to puncture the catheter shaft for urine sampling.
- There is a danger of needle stick injury and the catheter may be compromised.
- Do not handle with sharp instruments.
- This product is for single use only. Do not re-sterilize or reuse.

STORAGE

Do not store in extreme temperature or wet conditions. Contents considered sterile unless package is opened, damaged or wet.

This product was sterilized using a validated and controlled ETO sterilization process. Resterilization and reuse is prohibited and will cause contamination to the patient or the caregiver. The device should be disposed after use in accordance with the hospital policy for biological waste. Store in a cool dry place.

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