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GB

Please read the instruction manual carefully before using it!

The skin stapler which designed and produced by our factory is the special instrument for surgical skin suture. With novel structure and application security, our skin stapler has the advantages of quick speed of sewing up during the operation, little histological reaction, good wound matching and healing, small surgical scar, and no pain for removing, and so on. Our disposable skin stapler enables you to not only work with higher efficiency, but also avoid iatrogenic cross-infection effectively.

1. Main Structure

- a. The skin stapler is composed of suture instrument and stapler. Every skin stapler has 35pcs staplers inside. And the skin suture instrument is composed of the parts of body, handle, stapler box, stapler push device and stapler turn back device.
- b. The material for our stapler is 00Cr18Ni14Mo3 which is one medical used stainless steel. The material for the skin stapler body, handle & stapler box is ABS & plastic. And the stapler turn back device is made of 1Cr18Ni9Ti material.

2. Working Principle

For the aim of skin suture, after the handle was pushed by outside force, the staplers will be shaped after extrusion under the mechanical transmission of the each part.

3. Capability Index


- a. The suture force of the skin stapler should be no more than 55N
- b. The metal parts of the skin stapler have the good function of anti-eroded.
The skin stapler is sterilized with ethylene oxide.
- c. The skin stapler is sterilized with ethylene oxide. 14 days after the ETO sterilization, the leftover of the ethylene oxide should be less than 10 mg /kg.

4. Applicative area

Disposable skin stapler is used to suture incisions during surgical operation.

5. How-to-use

- a. Pull up the skin of both sides of the operative incision; clamp it with tissue tweezers so that both sides are aligned to each other.
- b. Make sure the bottom of the stapler is at the suitable angle to the middle line of the incision and the staple is vertical to the longitudinal axis of the incision, then put the stapler gently onto the incision till the fore-arrow levels with the middle line of the incision.
Caution: Don't put the stapler tightly to the incision; otherwise it will make removing staples difficult!
- c. Clasp the handle and don't loosen it until you hear a "click", then draw the stapler a little back and get ready for next manipulation.

6. Important Notice:

- a. The skin stapler must be operated by the professional person who has the medical certificates.
- b. The skin stapler is sterilized with ethylene oxide.
- c. The skin stapler could not be polluted when it was taken out of the package. Forbidden to use if the package broken.
- d. The product is limited to disposable use; please soak it in disinfectant according to Infection Control Criteria of your hospital before destroying/disposing.
- e. Disposable skin stapler is only used to suture incisions, if the distance of the skin incisions and the internal important histiocyte is less than 5mm, forbidden to use the skin stapler to do the skin suture.
- f. If the hypodermic part under incision is too thick, suggest to do the hypodermic suture first, then do the skin suture again.

- g. The suitable distance between one stapler to another should be between 5~10 mm. Please use the special skin remover to remove the staplers after 7~10 days. Or please refer to the doctor's prescribing.
- h. Please use the skin stapler within the period of validity.

7. Reserve and Transportation

The product should be kept in rooms with relative humidity $\leq 80\%$, with no erosive gases and good ventilation. Not to be stowed under other cargo, and handle with care.

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P. J. Dahlhausen & Co. GmbH