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Mölnlycke® Z-Flo™ Fluidized positioner

Z-Flo™ is a mark owned by EdiZONE, LLC of Alpine, Utah

Mölnlycke® Z-Flo™ Fluidized positioner

Mold the positioner to support anatomies or create boundaries

The positioner maintains its shape, securing desired patient position

en

Please read all of the instructions before use. Use as directed by appropriate clinician. Always follow your facility protocol for turning and positioning.

Mölnlycke® Z-Flo™ Fluidized positioner Neonatal/Pediatric

Positioner product for Neonatal/Pediatric in hospital population.

Product Description

The positioner consists of a fluidized medium composed of silicone oil and microspheres which is encapsulated in a polyurethane film exterior. The fluidized positioner can be molded to support multiple anatomical shapes and clinical needs.

Intended use

Mölnlycke® Z-Flo™ Fluidized positioner is designed to provide patient positioning, pressure redistribution and/or to offload anatomical sites, which may be used as part of a prophylactic therapy to help prevent pressure ulcers.

Indications for use

Neonatal/Pediatric

- Developmentally-appropriate infant positioning
- Patient positioning, turning and repositioning
- Off-loading of heel, sacrum, occiput and other bony prominences
- Contracture management

Lifetime

Product properties will deteriorate over time in use. The rate of deterioration is dependent on use environment and application. The product can be used with one patient until the product has deteriorated and does no longer fulfil its intended use.

Other information

If the product is used after expiry date product properties cannot be ensured. The product expiry date is printed on the packaging labels.

Instructions for use

GENERAL

The positioner can be molded to:

- achieve and maintain optimal position, cradling and creating boundaries for the patient.
- maximize surface contact area to minimize excessive forces on the skin.
- accommodate medical equipment or airway clearance.

The positioner can easily be re-molded to accommodate changed clinical needs. The positioners are safe to use in MRI, CT and X-ray. If X-ray image is taken through the fluidized positioner, wrinkles in film and cover should be avoided/minimized due to risk of image artefacts.

NEONATAL/INFANT

1. Remove the positioner from the transparent plastic bag.
- If required by the patient's clinical needs (e.g. in incubator), prewarm the positioner to desired temperature in a blanket warmer or equivalent.
3. Positioner cover or linen is required between the patient's skin and the fluidized positioner at all times during use.
4. Prior to placing the patient on the positioner, it can be pre-molded to accommodate certain anatomies or positions e.g. pre-molded to support the head or when placing patient in prone position.
5. Place the positioner appropriately under the patient.
6. Gently stabilize the patient on the positioner. Holding the patient with one hand, slowly mold the fluidized medium of the positioner into desired shape, supporting limbs and securing position of the patient.
7. If the product has deteriorated (e.g. due to time in use) or the size of the product no longer fulfils the clinical needs of the patient it should be discarded and replaced.

i

<https://www.molnlycke.com/zflo/neo>

Neonatal / Pediatric

PEDIATRIC

1. Remove the positioner from the transparent plastic bag.
2. Positioner cover or linen is required between the patient's skin and the positioner at all times during use.
3. Place the positioner appropriately under the patient.
4. Gently mold the positioner towards the patient/anatomical site to support and maintain the desired position, achieving pressure redistribution to minimize excessive forces on the skin.
5. If the product has deteriorated (e.g. due to time in use) or the size of the product no longer fulfils the clinical needs of the patient it should be discarded and replaced.

CLEANING AND MAINTENANCE INSTRUCTIONS:

- Use hospital approved disinfectant or household antimicrobial cleaner/disinfectant according to its label instructions. Avoid using cleaning agents containing bleach. Extended use of alcohol- based disinfectants may reduce the life of the product.
- If the fluidized medium feels non-homogeneous, the product can be restored by repetitive molding and shaping.

Warnings

- Patients under one year of age must be properly supervised and monitored in health care facilities when using positioners.
- Do not allow patients to lie face down in fluidized positioner unless closely monitored due to risk of suffocation.
- To avoid puncturing, be careful when handling sharp objects in proximity of the product.
- When handling the product, especially during molding, do not use fingernails and avoid wearing rings with raised or sharp edges.

Precautions

- The product should not be exposed to temperatures higher than 40°C (104° F).
- Do not autoclave or microwave.
- Replace promptly if punctured or damaged.
- If punctured, the positioner should be placed in a sealed plastic bag to avoid spillage of the content, which may increase risk of slips and falls.
- Closely monitor patient in prone position to secure airway clearance.

Special storage conditions/handling conditions

Mölnlycke Z-Flo Fluidized Positioners should be stored below 40°C (104°F). Long exposure to high temperatures may deteriorate the product prematurely.

Other info:

- Do not use on multiple patients (single patient use only).
- Product provided non-sterile. Do not sterilize.
- Not user repairable.
- Non-flammable according to SS-EN ISO 12952-1 & SS-EN ISO 12952-2.
- This product is non-toxic.
- Disposal should be handled according to local environmental procedures. The product can be incinerated.

If any serious incident has occurred in relation to the use of Mölnlycke® Z-Flo™ Fluidized positioner it should be reported to Mölnlycke Health Care and to your local competent authority.

Manufacturer

Mölnlycke Health Care AB
Gamlestadvägen 3C, Box 13080, SE-402 52 Göteborg, Sweden

The Mölnlycke trademark, name and logo are registered globally to one or more of the Mölnlycke Health Care Group of Companies

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