Mölnlycke®
Mepi™ Debripad

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Mepi™ Debripad

Microfiber pad for wound bed preparation

(en) INSTRUCTIONS FOR USE

Intended Purpose
Mepi™ Debripad is intended for the debridement of superficial, chronic and acute wounds and the skin surrounding the wound. It is to be used to remove devitalised tissue from the wound bed, and to absorb exudate and debris during debridement. It is suitable for dry, parchment-like, hyperkeratotic, seborrhoeic or oily healthy or damaged skin.

Product Description
Mepi™ Debripad is a polyester/polyamide microfiber pad designed for fast, safe and low pain debridement of wounds. Each pad is provided sterile and is for single use only.

Not made with natural rubber latex

Indications
Mepi™ Debripad is indicated for use on diabetic ulcers, arterial and venous ulcers, pressure ulcers, postoperative wounds healing by secondary intention, traumatic or surgical wounds (e.g. abrasions, incisions, lacerations), and, burns and scalds (3rd degree: following surgical debridement).

Contraindications
- Do not use on patients with a known sensitivity to any of the composition materials.
- Surgical implantation.
- Do not use on hard necrotic tissue or viable/granulating tissue.

Precautions or Warnings
- For use under the guidance of a healthcare professional.
- Single use only, do not re-use in whole or in part, as it may compromise sterility and/or the performance of the pad.
- Do not re-sterilise.
- Do not use if the package is damaged or if there are signs of contamination on the pad.
- Do not cut Mepi™ Debripad.
- Do not leave Mepi™ Debripad on the wound.
- Mepi™ Debripad must not be ingested and must be kept away from children and animals.
- Do not use if primary packaging is damaged or unintentionally opened before use.
- For disposal of the device please refer to local guidance on the disposal of clinical waste.

DIRECTIONS FOR USE

When using the Debridement pad standard practice should be followed, including the use of an aseptic technique, and sterile surgical gloves.

1. Remove any wound dressing.
2. Debridement Procedure:
 a. Soak Mepi™ Debripad in 20ml sterile saline solution or any other suitable wound cleansing solution.
 b. Insert hand into the mitt of the pad.
 c. Wipe the wound surface gently without pressure using the soft side of the pad that faces away from the palm of the hand.
 d. If required cleanse the peri wound skin using another thoroughly moistened Mepi™ Debripad.
3. Following use, standard treatment for the condition should be continued
4. Dispose of the used Mepi™ Debripad in line with local procedures for clinical waste.

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Date	14/11/2024 rev. C vj
Customer details	
Customer	Molnlycke
Product description	Mepi Debripad
Customer REF	300111, 300113
Finished Good	10017283
	N/A
Nr. of colours 1	BLACK Dieline

AMS Studio, Advanced Medical Solutions Ltd., Road One,
33 Premier Park, Winsford Industrial Estate, Winsford, Cheshire CW7 3RT UK
E-mail: ams.studio@admedsol.com Tel:+44 (0) 1606 86 3500

Supplier details	
AMS item no.	10017283 101/4/F7
Flat size	630mm (L) x 430mm (W) ± 1mm
Folded size	90mm x 70mm ± 2mm
Print Process	Offset litho
Raw Material ref	RMS006 - 60gsm
Front view when folded	Molylycke Mepi Debripad logo and title.

SUPPLIER
• Please acknowledge receipt of artwork as this legend will be the final specification and should match the RMS-006 Supplemental Supplier document.
• Check details against this the legend information should be correct and inline with the artwork.
• The supplier must check if this design can be printed on the substrate indicated.
• Before printing make sure that this legend matches the official order as well.
• Proceed to repro and proof for approval.
• Do not print until you receive an official order
• Do not destroy any previous plates until AMS confirm it is ok for you to do so. Any plates destroyed without this confirmation will be replaced at the supplier's cost.

Distributed by: Mölnlycke Health Care AB

The Mölnlycke trademark, name and logo are registered globally to one or more of the Mölnlycke Health Care group of companies

RESORBA Medical GmbH
Am Flachmoor 16, 90475 Nürnberg, Germany

CE 1434

2024-07 (101/4/F7)