

## Katétr AA63++ FOLYSIL® TIEMANN celosilikonový

### Výrobce

COLOPLAST A/S

### Distributor

COLOPLAST A/S

### CE třída

IIb (pravidlo 5)

### Sterilizace

Sterilní a jednorázové použití  
Etylénoxid (EO)

### Indikace

Katetrizace močového měchýře.  
Silikonové katétry Folysil se nesmějí ponechávat zavedené déle než 30 dní.

### Kontraindikace

Stejně jako u katetrizace močového měchýře a obecně veškeré známé alergické reakce na materiál použitý u výrobku (ZP).

*Návod k použití SH2115*

### Odkazy

CH / FR	Barva kroužku a ventilu	Balónek (ml)	Celková délka (cm)	Objednací čísla
12	bílá	10	40,2	AA6312
14	zelená	10	40,5	AA6314
16	oranžová	15	40,7	AA6316
18	červená	15	41	AA6318
20	žlutá	15	41,2	AA6320
22	fialová	15	41,5	AA6322
24	modrá	15	41,7	AA6324

### Popis

Tolerance délky: L <sup>+0</sup> /-25  
 Distální špička: coudé olivovité zakončení  
 Proximální konec: zúžený, kroužek a barevný ventil  
 Kanály: 2 (drenážní, pro nafukování)  
 Otvor: 1 uvnitř coudé špičky  
 Radioopakní hadička

**Materiály:** 100% silikon

**Obsahuje latex:** Ne

### Skladovatelnost

5 let

### Balení

**Uživatelské balení** 1 katétr v sáčku a sáčku typu peel pouch

**Objednávané balení** 5 uživatelských balení v karton. krabici 90x40x550 mm

Regulatorní záležitosti

Stéphane Bouché  
Regulatory Affairs Specialist  
Porges - a Coloplast division

Marketing  
François Guitteny  
International Brand Manager HCM  
Porges - a Coloplast Division




## AA63++ FOLYSIL® TIEMANN catheter

### All silicone

**Manufacturer**

COLOPLAST A/S

**Distributor**

COLOPLAST A/S

**CE class**

IIb (Rule 5)

**Sterilisation**

 Sterile and single use  
 Ethylene oxide (EO)

**Indications**

 Urethral urinary catheterisation.  
 Folsil silicone catheters must not be left in place for more than 30 days.

**Contraindications**

Same as urethral urinary catheterisation, and generally, known allergic reactions due to the device material).

*Instructions for use SH2115*
**References**

CH / FR	Ring and valve colour	Balloon (ml)	Total Length (cm)	Order references
12	White	10	40.2	AA6312
14	Green	10	40.5	AA6314
16	Orange	15	40.7	AA6316
18	Red	15	41	AA6318
20	Yellow	15	41.2	AA6320
22	Purple	15	41.5	AA6322
24	Blue	15	41.7	AA6324

**Description**

 Length tolerance: L <sup>+0</sup> / <sub>-25</sub>  
 Distal tip: coudé olive tip  
 Proximal end: funnel, ring and coloured valve  
 Channels: 2 (drainage, inflation)  
 Eye: 1 inside the coudé tip.  
 Radiopaque line

**Materials:** All silicone

**Containing latex:** No

**Shelf life**

5 years

**Packaging**

<b>Use unit</b>	1 catheter in a pouch and peel pouch
<b>Order unit</b>	5 use units in a carton box (mm) 90x40x550

Marketing

 François Guitteny  
 International Brand Manager HCM  
 Porgès - a Coloplast Division



Regulatory Affairs

 Stéphane Bouché  
 Regulatory Affairs Specialist  
 Porgès - a Coloplast division

